



AMG
EQUIPMENT



SERVICE DESIGNED WITH THE BUYER IN MIND

WE PROVIDE:

- | | |
|------------|-------------|
| Case Carts | Pads |
| Tables | Supports |
| IV Poles | Positioners |
| Hampers | Mattresses |
| Warmers | Leg Holders |
| Cushions | ...and more |



AMG is a healthcare equipment distribution company servicing hospitals and healthcare providers nationwide.

Our objective is to provide exceptional service while providing healthcare buyers access to a comprehensive portfolio of products across several trusted manufacturers. Operated and owned by former hospital purchasing and supply chain personnel, we allow buyers to fulfill several projects across multiple departments with one call or e-mail.

Our story is simple.

AMG was founded by former hospital employees with 20+ years in healthcare supply chain experience with a mission

to help buyers navigate the crowded market of suppliers. Having been in your shoes before, we know how to exceed your expectations and minimize your day-to-day challenges.

Whether it is a single routine equipment purchase, or a new construction project, AMG is focused on making the request for quote process simple and seamless. E-mail our customer service department with your need(s) and/or simply attach your equipment list and AMG will tap into our comprehensive portfolio of manufacturers and models to promptly return a quote with competitive pricing.

With AMG's ability to look across multiple manufacturers covering multiple departments, we save you time and money. Contact us today and find out for yourself.

Index

Case Carts.....	4	Armboards.....	20
IV Poles.....	6	Leg Support Pads.....	22
Hampers.....	8	Gel Pads & Positioners.....	24
Tables & Warmers.....	10	Safety Straps & Positioners.....	26
Waterloo Carts.....	12	Mattresses.....	28
Cushions.....	14	Leg Holders & Accessories.....	30
Table Pads.....	16	Siderail Clamps & Table Extenders.....	34
Arm Support Pads.....	18	Refurbished Items.....	36

Case Carts

CDS-233

Cabinet Dimensions

23"W x 24¼"D x 47"H

Usable Interior Space

20¼"W x 23½"D x 44½"H

Overall Dimensions

28½"W x 28½"D x 56"H

Shown with optional shelves

CDS-245

Cabinet Dimensions

40¼"W x 24¼"D x 47"H

Usable Interior Space

37½"W x 23½"D x 44½"H

Overall Dimensions

45½"W x 28½"D x 56"H

Shown with optional shelves and cases

RCC-233

- Single door
- re•volve™ 6" diameter, maintenance-free, high-temp, engineered thermoplastic casters
- Two swivel and two with direction lock

Overall outside dimensions

29"W x 27½"D x 59"H

Usable Interior Space

20"W x 23½"D x 44¼"H

Shelving Configurations:

RCC-233-SS

Solid shelf, for RCC-233 and RCC-256-MS.

RCC-233-SSRO

Roll-out solid shelf, with (2) position spring loaded shelf lock. Shelf locks in storage and fully extended position for RCC-233 and RCC-256-MS.

RCC-233-WS

Wire shelf, for RCC-233 and RCC-256-MS.

RCC-233-WSRO

Roll-out wire shelf, with (2) position spring loaded shelf lock. Shelf locks in storage and fully extended position for RCC-233 and RCC-256-MS.

RCC-256-MS

- Single door
- re•volve™ 6" diameter, maintenance-free, high-temp, engineered thermoplastic casters
- Two swivel and two with direction lock

Overall outside dimensions

29"W x 27½"D x 42"H

Usable Interior Space

20"W x 23½"D x 27¾"H

Shelving Configurations:

RCC-233-SS

Solid shelf, for RCC-256-MS and RCC-233.

RCC-233-SSRO

Roll-out solid shelf, with (2) position spring loaded shelf lock. Shelf locks in storage and fully extended position for RCC-256-MS and RCC-233.

RCC-233-WS

Wire shelf, for RCC-256-MS and RCC-233.

RCC-233-WSRO

Roll-out wire shelf, with (2) position spring loaded shelf lock. Shelf locks in storage and fully extended position, for RCC-256-MS and RCC-233.

CCC1-17

Maxi Case Cart comes with (1) full length heavy duty stainless steel wire roller shelf. The CCC1-17 is recommended for procedures requiring a larger volume or size of instruments and/or supplies, such as neurosurgery and orthopedic surgery. **Body** – 18 gauge, **Door(s)** – 18 gauge, **Top** – 16 gauge, **Chassis** – 11 gauge, **Bumper channel** – 14 gauge.

CCC2-17

Multi-Purpose Case Carts are the right size for most routine surgical procedures, whether or not your facility utilizes a container system. **Body** – 18 gauge, **Door(s)** – 18 gauge, **Top** – 16 gauge, **Chassis** – 11 gauge, **Bumper channel** – 14 gauge.

CCC2E-17

Multi-Purpose Case Carts are the right size for most routine surgical procedures, whether or not your facility utilizes a container system. **Body** – 18 gauge, **Door(s)** – 18 gauge, **Top** – 16 gauge, **Chassis** – 11 gauge, **Bumper channel** – 14 gauge.

CCC3-17

Mini Case Cart is designed for labor and delivery, outpatient surgery or whenever a smaller cart is preferred to transfer instruments and supplies from one area to another. **Body** – 18 gauge, **Door(s)** – 18 gauge, **Top** – 16 gauge, **Chassis** – 11 gauge, **Bumper channel** – 14 gauge.

CCC4-17

Space-Saver Case Carts are designed to give added capacity by utilizing the cube while minimizing the floor space required by larger carts. **Body** – 18 gauge, **Door(s)** – 18 gauge, **Top** – 16 gauge, **Chassis** – 11 gauge, **Bumper channel** – 14 gauge. Removable pull-out type stainless steel wire shelves are provided on closed carts. Shelves slide on adjustable angle sliders, suspended from keyhole punched standards welded to the cart interior. Open carts offer either a solid stainless or wire shelf, adjustable on 4" centers.

CCC5-17

Ultra Space-Saver Case Carts offer greater capacity using the same floor space as the Multi-Purpose Cart. It's useful for routine or extensive procedures requiring increased volume/size of instruments and/or supplies. **Body** – 18 gauge, **Door(s)** – 18 gauge, **Top** – 16 gauge, **Chassis** – 11 gauge, **Bumper channel** – 14 gauge. Removable pull-out type stainless steel wire shelves are provided on closed carts. Shelves slide on adjustable angle sliders, suspended from keyhole punched standards welded to the cart interior. Open carts offer either a solid stainless or wire shelf, adjustable on 4" centers.




CDS-233



CDS-245



RCC-233



RCC-256-MS



CCC1-17



CCC2-17



CCC2E-17



CCC3-17



CCC4-17



CCC5-17

Roller Board Wall Mount

PT-5800

Store roller boards on the wall
Durable stainless steel mounting brackets keep transfer devices up off the floor.
Instructions for mounting, placement, and hardware are included.

Patient Roller Board Transfer Devices

Transfers patients from one surface to another with minimum effort.
Reduces the risk of back injuries.
Aluminum rollers allow easy glide.
Heavy-duty vinyl cover – easy to clean.
Replacement vinyl covers available.

PT-5920

20"x 31" Heavy duty construction
Ideal for Bariatric patients.

Surgical Scrub Sinks

All stainless steel construction
Manual or automatic operation
Adjustable thermostatic water mixing valve
Single, double, and triple sink basin.
Optional Features Available.

IV Poles

IVP-4200

S.S. IV Pole w/ 2 hook, extends 53"- 96"
5 leg powder coated base.

IVP-4400

S.S. IV Pole w/ 4 hook, extends 53"- 96"
5 leg powder coated base.

IVP-4600

S.S. IV Pole w/ 6 hook, extends 53"- 96"
5 leg powder coated base.

IVP-6050

Chrome heavy-duty IV Pole, 5 leg base,
Tru-loc knob, 4 hook top, extends 56"-100"

IVP-6250

S.S. 4 leg IV Pole w/ 2 hook top,
Tru-loc knob, extends 51"-94"H

IVP-6450

S.S. 4 leg IV Pole w/ 4 hook top
Tru-loc knob, extends 51"-94"H

IV Poles Accessories

Economy IV Poles

IVP-1200

Chrome IV Pole, 4 leg, 2 hook 47"- 85"

IVP-1400

Chrome IV Pole, 4 leg, 4 hook 47"- 85"

Versa-Pole

IVP-8800

Portable I.V. pole w/clamp secures where you need it. Works on stretchers, wheelchairs, etc.
S.S., 2 hook, 22"- 36" high

Super Pole Unit

IVP-8900

Stainless Steel telescoping pole,
2 hook top, 33"- 60",
Securely mounts on beds, gurneys,
and medical carts.

IV Pole Accessories

IVP-1050

S.S. Infusion pump tray/utility tray for I.V. Pole.

IVP-1090

Infusion pump support frame with
2 universal pole clamps for mounting.

IVP-8910

Patient Pal Support Wheel
Holds 4 utility trays.

IVP-8915

Utility Tray Triangle
Durable plastic, 4 per package.

IVP-8950

Standard C-Clamp for IV Poles.

EZ Power


IVP-9080

6 receptacle outlet hub
Securely mounts on any IV Pole.
15' power cord, 120v. 15 amp.

IVP-9085

6 receptacle outlet hub
Securely mounts on any IV Pole.
15' power cord, 120v. 20 amp.

IV Poles



IVP-1200



IVP-4200



IVP-6050



IVP-6250



IVP-7140



IVP-7240



IVP-8800



IVP-8900

Roller Board Wall Mount



PT-5825
Shown on PT-5800

PT-5800

Patient Roller Board Transfer Devices



PT-5920

Surgical Scrub Sinks



EZ Power



IVP-9080

IV Poles

Series 1752

P-1572-4

4 hook

P-1572-CV

S.S. IV Pole w/2hook extends 53
5 leg powder coated.

Series 1576

P-1576-4

4 hook

P-1576-CV

6 hook

P-1080-6

Designed specifically for ambulatory
infusion therapy applications.

P-1083-6

Designed specifically for ambulatory
infusion therapy applications and
made from CuVerro™ bactericidal
copper alloy.

Series 1076

P-1076-6

6 hook

P-1080-6

Designed specifically for ambulatory
infusion therapy applications.

P-1083-6

Designed specifically for ambulatory
infusion therapy applications and
made from CuVerro™ bactericidal
copper alloy.

SERIES 1076

P-1076-6

6 hook

Hampers

Newly designed lidded hampers
feature the Quiet-Close™ lid dampener
for improved safety and added control
of airborne bacteria and odors. Linen
hampers are constructed of heavy
gauge steel and feature a triangular
shape for efficient use of space. Double
ball bearing 3" swivel casters make
for easy mobility. Strong, heavy gauge
shelf keeps bag from touching floor
and simplifies bag removal.

P-120

32"H with a 59½" perimeter

P-121

35"H with a 80¾" perimeter

P-1120-SS

32"H with a 59½" perimeter

P-1121-SS

35"H with a 80¾" perimeter

Features/Specs

- Powder coated steel
- Accommodates 18" bags
- Open Linen Hamper

P-120-L

32"H with a 59½" perimeter

P-121-L

35"H with a 80¾" perimeter

P-1120-L-SS

32"H with a 59½" perimeter

P-1121-L-SS


35"H with a 80¾" perimeter

Features/Specs


- Stainless steel
- Includes rolling bumpers
- Foot operated
- Accommodates 18" bags
- Quiet-Close™ Lid



P-120
P-121
P-1120-SS
P-1121-SS



P-120-L
P-121-L
P-1120-L-SS
P-1121-L-SS



* approximately the size of a basketball



IV Poles




P-1572 SERIES




P-1576 SERIES



P-1080-6



P-1083-6



P-1076 SERIES

Foot pedal height adjustment

19" Caster Center Line
21½" Footprint

Tables

CDS-1100

Main Table Top Surface: 72"W x 30"D
Height Adjustable: 32" to 43 1/4"H

Mid-Size Tables

WITH A SHELF

Dimensions:

SG-82-SS 20"W x 16"D x 34"H

SG-84-SS 30"W x 16"D x 34"H

SG-86-SS 33"W x 18"D x 34"H

SG-88-SS 36"W x 20"D x 34"H

WITHOUT A SHELF

Dimensions:

SG-83-SS 20"W x 16"D x 34"H

SG-85-SS 30"W x 16"D x 34"H

SG-87-SS 33"W x 18"D x 34"H

SG-89-SS 36"W x 20"D x 34"H

Large Tables

WITH A SHELF

Dimensions:

SG-90-SS 48"W x 20"D x 34"H

SG-92-SS 48"W x 24"D x 34"H

SG-94-SS 60"W x 24"D x 34"H

SG-96-SS 72"W x 24"D x 34"H

SG-98-SS 36"W x 24"D x 34"H

WITHOUT A SHELF

Dimensions:

SG-91-SS 48"W x 20"D x 34"H

SG-93-SS 48"W x 24"D x 34"H

SG-95-SS 60"W x 24"D x 34"H

SG-97-SS 72"W x 24"D x 34"H

SG-99-SS 36"W x 24"D x 34"H

Instrument Tables

7830SS

Instrument Table with shelf offers complete 18 gauge, 304 stainless construction. Lightweight, yet extremely durable. Sanitary closed top hem prevents bacterial build-up. Sealed reinforcing channels provide structural rigidity. 1 1/4" diameter legs are fully welded for strength and stability. There is a 16" clearance between shelves on the Howard Models.

3026SS

Sawyer Instrument Table with H-Brace features a new and improved design offering complete 18 gauge, 304 stainless construction. Lightweight, yet extremely durable. Sanitary closed top hem prevents bacterial build-up. Sealed reinforcing channels provide structural rigidity. 1 1/4" diameter legs are fully welded for strength and stability. The Sawyer Table has an H-brace in lieu of a second shelf. H-braces are fully welded around the legs with no crevices for dirt to collect. All tables travel on three-inch casters; shipping weight 40 lbs.

Warmers

P-2140

Blanket & Fluid Warmers

Deluxe blanket warmers feature 4 digit LED display, window doors, control lock programming and auto time programming.

P-2140

OVERALL SIZE
• 25 1/4"W x 35"D x 72 3/4"H
COMPARTMENT SIZE
• 7.7 cubic ft. blanket compartment
• 18 7/8"W x 27"D x 26 1/8" H
• 7.7 cubic ft. fluid compartment
• 18 7/8"W x 27"D x 26 1/8" H

- Top compartment for blankets, lower compartment for fluids
- Fluid compartment has three epoxy-coated baskets each with a maximum capacity of 24 liters. Total capacity is 72 liters
- Separate temperature controls for each compartment
 - Temp Range Blankets: 98°F - 200°F
 - Temp Range Fluids: 98°F - 150°F
- Perimeter bumper included

P-2148

Blanket & Fluid Warmer

Deluxe blanket & fluid warmers feature 4 digit LED display, window doors, control lock programming and auto time programming.

P-2148

OVERALL SIZE
• 31"W x 21 1/4"D x 66 1/2"H
COMPARTMENT SIZE
• 28 liter fluid compartment
• 26 1/2"W x 17 1/8"D x 12 1/2" H
• 8.7 cu. ft. blanket compartment
• 26 1/2"W x 17 1/8"D x 32 1/2" H

- Top compartment for fluids, lower compartment for blankets
- Fluid compartment has one easy access epoxy-coated basket with a maximum capacity of 18 1-liter bottles or equivalent
- Separate temperature controls for each compartment
 - Temp Range Blankets: 98°F - 200°F
 - Temp Range Fluids: 98°F - 150°F
- 18" depth perfect for built-in OR applications

P-2205

P-2205

OVERALL SIZE
• 18 1/2"W x 19 7/16"D x 16 1/8"H
COMPARTMENT SIZE
• 1.5 cubic ft.
• 14 1/2"W x 15 1/2"D x 9 7/8"H

- Unit can be stacked
- Temp Range: 90°F - 160°F

P-2240

P-2240

OVERALL SIZE
• 24"W x 27"D x 34 13/16"H
COMPARTMENT SIZE
• 7.5 cubic ft.
• 21"W x 20 15/16"D x 28 1/4"H
• (13 3/4"H above and below the center shelf)

- Unit can be stacked
- Temp Range: 90°F - 160°F
- Heated center shelf

Tables

OPTIONS

Accessory

- Work Lamp CDS-1020**
 - Work Lamp arm has a 42" reach
 - Features a 5" viewing lens
 - 75% magnification
- Outlet Strip CDS-1015**
 - Hospital grade
 - Six outlets
 - 15' foot cord
- Keyboard Tray CDS-1032**
 - Measures 22 1/4"W x 10"D x 2 1/4"H
 - Tray extends 10 1/2" (retractable)
 - Built-in mouse tray comes standard
- Computer Enclosure CDS-1031**
 - Open back for air flow and easy access to computer cords.
 - Enclosure can be customized to accommodate larger computer towers (max. load 30 lbs/13.60 kg)
- Extension Arm CDS-1030**
 - Supports LCD monitors (max. load 25 lbs / 11.3 kg)
 - For optimal viewing, extension arm features multi-directional adjustment
 - Gas-strut and ball-joint system for smooth, stable operation
- Storage Cart CDS-1010**
 - 4 drawer tote box cart for added storage fits under the workstation
 - Totes are constructed of impact-resistant, high-density polyethylene and measure 17 3/4"L x 14 1/2"W x 5"D

CDS-1100



7830SS



3026SS



SG-82-SS
SG-84-SS
SG-86-SS
SG-88-SS



SG-83-SS
SG-85-SS
SG-87-SS
SG-89-SS



SG-90-SS
SG-92-SS
SG-94-SS
SG-96-SS
SG-98-SS



SG-91-SS
SG-93-SS
SG-95-SS
SG-97-SS
SG-99-SS

Warmers



P-2140



P-2148



P-2205



P-2240



DURABLE finishes.

All Waterloo carts are coated with a highly durable powder paint that is designed to withstand the demands of the medical environment. The powder paint we use is designed to stand up to standard cleaning solvents used in central sterile supply for surface cleaning.

With 16 colors to choose from and the ability to mix and match drawer colors on a single cart you will be able to create a unique cart for your particular application and/or to match a specific decor.

Color Options:



Recommended Accessories



Tilt Bins
TB-3
Clear plastic tilt bins with white exterior. Available in 3, 4, 5, or 6 bin units. Mounts to base rail system (BRS-2, BRS-3 OR BRS-JR)
TB-3
TB-4
TB-5
TB-6



Small Clear Plastic Organizer
ORG-1S
Includes clear slide-in dividers. Mounts to accessory railing.
Size:
9"W x 5"D x 3"H



Horizontal Rail System
BRS-2
For hanging accessories and mounting shelves on Unicarts. Includes two vertical and three horizontal rails. For Steel or Aluminum Unicarts purchase BRS-2. Med Jr. carts purchase BRS-JR.



Tape Dispenser
TD-1
White tape dispenser with suction cups to attach to the top of a cart.



Large Waste Basket
BSK-2
Large (20 quart) plastic waste basket. Requires accessory rail.
Size:
13"W x 8½"D x 15"H



Extra-Large Wire Basket
WSB-1
Mounts to accessory railing.
Size:
15½"W x 9½"D x 5½"H

Emergency Carts

Rapid and organized response to a cardiac emergency is essential. The Waterloo aluminum cart is about 30% lighter than a standard steel cart, offering faster response times. The included steering caster allows the cart to take tight corners further improving maneuverability.

Work surface and drawer organization are vital during a Code. Our dual slide-out shelves virtually double your cart top work surface, allowing you to move equipment/supplies to the cart top and within easy reach.



UXRLA-333669-RED 6-Drawer X-Tall
Dimensions:
29"w x 24½"d x 45"h



UMRLA-33369-RED 5-Drawer Medium
Dimensions:
29"w x 24½"d x 39"h



USBLA-3369-ELB 4-Drawer Short
Dimensions:
29"w x 24½"d x 36"h

Anesthesia Carts

Economical steel carts offer a robust work platform that can be configured with accessories to address your anesthesia cart needs. Shares many of the same great features as our aluminum cart but at a more economical initial price. Other features include:

- Smooth ball bearing slides with auto-close mechanism to pull the drawer fully closed.
- Advanced nylon reinforced polycarbonate bumper provides a wider wheel base and greater cart stability and protection. Single piece construction resists scuffs and scratches.
- 5" Premium "Easy roll" sealed bearing casters with precision bearing swivels. Thread guards minimize sutures and other materials from collecting around axles and keep carts rolling smoothly.



USGKU-33339-TLG 5-Drawer Short
Dimensions:
29"w x 24½"d x 36"h



UTGKU-33369-DKB 5-Drawer Tall
Dimensions:
29"w x 24½"d x 42"h



UXGSU-333669-ELB 6-Drawer X-Tall
Dimensions:
29"w x 24½"d x 45"h


Surgical Table Cushion Sets

- Standard Series
- Pro-Tek Series
- Integra-Gel Series
- Bariatric Comfort Series



Pad Color Options:


Please specify choice of color when ordering. Unless specified, black vinyl will be supplied.



Standard Series

Recommended for low risk patients, short duration procedures. 2" thick overall. Constructed of medium firm high density foam, covered with 2 way stretch vinyl.


- 2-Way Stretch Vinyl
- Medium Firm High Density Foam



Pro-Tek Series

Designed for Pressure Management, high risk patients, medium to long duration procedures. 2 1/2" thick overall. Constructed of medium firm high density comfort foam, a top layer of anti-shear Pro-tek Memory foam, covered with 4-way stretch ultra soft vinyl.


- 4-Way Stretch Vinyl
- Pro-tek Memory Foam
- Medium Firm High Density Foam



Integra-Gel Series

Designed with ultimate patient comfort in mind. Recommended for high risk and heavy weight patients. 3" thick overall. Constructed of medium firm high density comfort foam, a layer of anti-shear Pro-Tek Memory foam, and a top layer of visco polymer Gel, covered with 4-way stretch ultra soft vinyl.

- 4-Way Stretch Vinyl
- Visco Polymer Gel
- Pro-tek Memory Foam
- Medium Firm High Density Foam



Bariatric Comfort Series


Designed with the comfort of bariatric patients in mind. Recommended for low - high risk bariatric patients. 4" thick overall. Constructed of high density firm foam, Medium Firm High Density foam, and a layer of anti-shear Pro-tek Memory foam covered with 4-way stretch ultra soft vinyl.

- 4-Way Stretch Vinyl
- Pro-tek Memory Foam
- Medium Firm High Density Foam
- High Density Foam

MAKE	MODEL	Standard Series	Pro-Tek Series	Integra Gel Series	Bariatric Comfort Series
AMSCO/ STERIS	2080M - 2080L - 2080RC - 1080	SC-3220	PC-3220	IGC-3220	BCS-3220
	3080RC - 3080RL - 3080SP - 3085SP	SC-3320	PC-3320	IGC-3320	BCS-3320
	4085/CMAX	SC-3380	PC-3380	IGC-3380	BCS-3380
	5085	SC-3390	PC-3390	IGC-3390	BCS-3390
	5085 SRT	SC-3395	PC-3395	IGC-3395	BCS-3395
	ORTHOGRAPHIC II (2-pc set) - (Optional Pads available)	SC-3540	PC-3540	IGC-3540	BCS-3540
	MAQUET	1132 Alpha Star (4-pc. Set)	SC-5380	PC-5380	IGC-5380
MIDMARK- CHICK	7100	SC-4550	PC-4550	IGC-4550	BCS-4550
	7300	SC-4560	PC-4560	IGC-4560	BCS-4560
SHAMPAINE- AFFILIATED- MDT	1501 - 1502 - 1503 - 1505	SC-1320	PC-1320	IGC-1320	BCS-1320
	RADI-OP 1900 - 1900E - 1900RC	SC-1520	PC-1520	IGC-1520	BCS-1520
	RADI-OP 4900E - 4900B - 5100	SC-1620	PC-1620	IGC-1620	BCS-1620
SKYTRON	ELITE 3100	SC-7030	PC-7030	IGC-7030	BCS-7030
	ELITE 3500	SC-7035	PC-7035	IGC-7035	BCS-7035
	3600B	SC-7040	PC-7040	IGC-7040	BCS-7040
	5001 - 6001 - 6001K - 6001KR - 6001KRB - 6500 - 6500B - 6500HD (3-pc. Set)	SC-7020	PC-7020	IGC-7020	BCS-7020
	6600	SC-7230	PC-7230	IGC-7230	BCS-7230
	6700B - 6701	SC-7280	PC-7280	IGC-7280	BCS-7280

*Replacement cushion sets are available for all makes and models of surgical tables, orthopedic tables, and other OR equipment.

- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial
- Latex Free



Waterfall Edge Construction



Pro-tek Memory Foam

Replacement Table Pads

Pad Color Options:

Please specify choice of color when ordering. Unless specified, black vinyl will be supplied.



Replacement OSI Table Pads

Recover existing pads except for CPS-9340, PC-6817, and MP-1013.

PRT-9100
Table Hip Pad Assembly - Large Patient

PRT-9130
Table Hip Pad Assembly - Small Patient

PRT-9150
Table Hip Pad Assembly - Regular Patient

PRT-9180-R
Thigh Pad - Right

PRT-9180-L
Thigh Pad - Left

PRT-9240
Chest Pad

PRT-9245
Dual Chest Pad

CPS-9340
Fusion Series Pad w/Hard Core - 6"W x 22"L x 1"T
2" Radius Corners - 16"L Loop Velcro
Sold individually

MP-1013
OSI Jackson sling with dual slot plastic hardware

PC-6817
Pro-tek Series replacement table pad for OSI Jackson 5803
21"W x 82"L x 3"T



PRT-9100



PRT-9130



PRT-9150



PRT-9180



PRT-9240



PRT-9245



CPS-9340



MP-1013



Arm Support Pads

Pad Color Options:

Please specify choice of color when ordering. Unless specified, black vinyl will be supplied.



Armboard Pads

(See table below)

1", 2", 3", or 4" pads available.

Flame retardant, Anti-microbial, and fluid resistant. Standard vinyl For 24" or 26" standard surgical table armboard. Wide Armboards and Pediatric Armboard pads available.

Arm & Hand Table Pads

Standard Style or Pro-Tek™

HT-1201-1

1" Standard Pad for HT-1201 Arm & Hand Table

HT-1601-1

1" Standard Pad for HT-1601 Arm & Hand Table

HT-1801-1

1" Standard Pad for HT-1801 Arm & Hand Table

HT-5201-2PM

2" Pressure Management Pad for HT-5201 'T' or Winged Arm & Hand Table

HT-6601-2PM

2" Pressure Management Pad for HT-6601 Rectangle Arm & Hand Table

HT-6901-2PM

2" Pressure Management Pad for HT-6901 Hourglass Arm & Hand table

HT-8201-2PM

2" Pressure Management Pad for HT-8201 Arm & Hand Table

Incline Armboard Pads

For all standard 26" armboards. Flame retardant. Anti-microbial and fluid resistant.

AS-4001-CT2

Incline from 2" to 6" thick

AS-4001-CT3

Incline from 3" to 7" thick

Vinyl Pads for Armguards

Memory foam pads covered with 4-way stretch Fabric. Easily slips over armguard. Disposable pads available (See pg 38)

AS-8020

Vinyl Soft Pad provides patient protection and comfort

AS-8040

Double-Care Pad provides protection for both patient and caregiver. Memory foam on both sides of armguard.

AS-8520

Vinyl Soft Pad provides patient protection and comfort. For Bariatric Armguard.

AS-8540

Double-Care Pad provides protection for both patient and caregiver. Memory foam on both sides of armguard. For Bariatric Armguard.

Multi-Purpose Armboard Pads

Pressure management pads – 2" thick, 4-way stretch fabric, Flame retardant, Anti-microbial, and fluid resistant. For multi-purpose surgical table armboards.

AS-6001-2PM

For 6" x 20" Multi-purpose Armboard

AS-6001-2CT Contoured

For 6" x 20" Multi-purpose Armboard

GP-2160

For 6" x 20" Multipurpose Armboard contoured pad w/gel overlay

Custom Pads

Custom replacement pads are available for all styles of arm supports. Specify foam style, fabric, thickness, and attachment device requirements when ordering.



AS-4001-2PM

Armboard Pads



AS-1201-2PM



AS-4001-2PM



AS-4111-2PM

Arm & Hand Table Pads



HT-5201-2PM



HT-6901-2PM



HT-1601-1

Armguard Pads



AS-8020



AS-8040

Incline Armboard Pads



AS-4001-CT2

Multi-purpose Armboard Pads



AS-6001-2CT



GP-2160



AS-6001-2PM

Custom Pads



Custom Armboard Pad



Custom Wide Armboard Pad



Custom Arm & Hand Table Pad

Size	Pad	Standard	Pro-Tek PM	Integra Gel
6" x 24"	1"	AS-2001-1	N/A	N/A
6" x 24"	2"	AS-2001-2	AS-2001-2PM	AS-2001-2ICG
6" x 26"	1"	AS-4001-1	N/A	N/A
6" x 26"	2"	AS-4001-2	AS-4001-2PM	AS-4001-2ICG
6" x 26"	2.5"	N/A	AS-4001-2.5PM	N/A
6" x 26"	3"	AS-4001-3	AS-4001-3PM	AS-4001-3ICG
6" x 26"	4"	AS-4001-4	AS-4001-4PM	AS-4001-4ICG
11" x 26"	2"	AS-4111-2	AS-4111-2PM	N/A
6" x 20" Pedi	2"	N/A	AS-1201-2PM	N/A

Standard Armboards

Phenolic operating surface. Rotates 180°, one-finger operation. 21 locking increments. Standard 24" or 26" length. Easy to attach to all North American surgical tables. Pressure Management and standard pads available in 1", 2", 3", or 4" (See pg 12)

AS-2001

6" x 24" Standard Armboard w/o pad

AS-4001

6" x 26" Standard Armboard w/o pad

Standard Armboard with Pad

AS-2000

24" Standard Armboard with 2" pad

AS-4000

26" Standard Armboard with 2" pad

AS-4000-1

26" Standard Armboard with 1" pad

AS-4000-3

26" Standard Armboard with 3" pad

AS-4000-4

26" Standard Armboard with 4" pad

Clamping Armboard

AS-4051

A spring loaded rail-lock clamping mechanism, allows for secure positioning on siderail during procedures requiring the table back to be raised. 6" x 26" Phenolic operating surface. Rotates 180°; 21 locking increments. Pads available in 1", 2", 3", or 4" (See pg 12)

Wide Armboard

AS-4110

11" x 26" phenolic operating surface. Easy one-finger positioning; 21 locking increments. Wide surface provides for ultimate patient comfort. Comes w/2" Pressure Management pad. Mounts to any North American surgical table siderail.

Duralux™ Armboard

AS-4501

Duralux™ operating surface, Standard 26" length. Rotates 180°; 21 locking increments. Easy to attach to all North American surgical tables. Pressure Management and standard pads available in 1", 2", 3", or 4" (See pg 12)

Quick Release Armboard

AS-4950

100% Radiolucent Carbon Fiber. 6" x 26" radiolucent surface. Quick release mounting device. Non-locking swivel gear allows for quick and easy 180° positioning. Comes w/2" Pressure Management pad. Mounts to any North American surgical table siderail.

Shoulder Mount Armboard

AS-4970

100% Radiolucent Carbon Fiber. 6" x 26" radiolucent surface. Slides under table top cushion. No mounting device, no metal parts. Non-locking swivel gear allows for quick and easy 180° positioning. Comes w/2" Pressure Management pad.

Armboard Holder

AS-5800

Wall Mount Holder for Standard Armboards. 6.5"W x 4"D x 15"H. Heavy duty anodized aluminum. Mounting hardware included.

Multi-Position Armboard

AS-6000

6"W x 20"L Phenolic operating surface. Permits complete range of positioning. Stainless steel support bar. Comes w/2" Pressure Management pad. Provides arm support for Lateral, Prone, Neuro, Arm Chair, and Park Bench positioning. Requires clark socket or super clamp for mounting to siderail of surgical table. For contoured armboard pad

Kross Arm Support

AS-7500

Perfect for arm chair, prone, lateral, and neuro procedure positioning. SS arm cradle measures 6" x 20". Disposable Kross Cloth sleeve makes this a sanitary choice, one roll included. Requires one clark socket for mounting (Kross Cloth available, See pg 38)

Surgical Table Armguard

AS-8000

Radiolucent. 8" High x 16" Long. Protects patient extremities. Easily slips under table cushion. Prevents nerve damage to patients. (Armguard Pads available - see pg 12, 38)

AS-8500

For Bariatric Patients. Radiolucent. 10.5"H x 16"L. Protects patient extremities. Easily slips under table cushion. Prevents nerve damage to patients. (Armguard Pads available - see pg 12)

Pediatric Armboard

AS-1200


6" x 20" Phenolic operating surface. Provides close efficient access for staff when working with small patients. Comes with 2" thick pad. Rotates 180°; 21 locking increments. Effortlessly mounts to siderail of any North American surgical table.

Standard Armboards



AS-4001
Pad not shown

Standard Armboard with pad



AS-4000

Clamping Armboard



AS-4051

Wide Armboard



AS-4110

Duralux™ Armboard



AS-4501

Quick Release Armboard



AS-4950

Shoulder Mount Armboard



AS-4970

Armboard Holder



AS-5800

Multi-Position Armboard



AS-6000-CT

AS-6000

Kross Arm Support



AS-7500

Surgical Table Armguard



AS-8000

Pediatric Armboard



AS-1200

Pad Color Options:

Please specify choice of color when ordering. Unless specified, black vinyl will be supplied.



Replacement Legholder Pads

Memory foam pads covered with 4-way stretch for comfort. Regular and Clamshell style available. Most styles listed are in stock. For all other makes and models of legholders, please call customer service.

ALS-6020
Altima™ Legholder Pads
Limited Stock. Old Style.

ALS-6025
(New) Altima™ Legholder - Clamshell Style Pads

ALS-6030
Altima™ Legholder Pads - Clamshell style
Limited Stock. Old Style.

LS-6220
Allen® black boots

LS-6230
Allen® black boots - Clamshell style

LS-9010
Ultrafin - Clamshell style

LS-9020
Yellofin™ boots - Clamshell style

LS-9020-D
Yellofin™ boots - Clamshell style
Disposable convoluted foam w/velcro

LS-9025
Yellow Banana boot



ALS-6025



LS-9501-2



Knee Crutch Pads

LS-9028
Skytron Levitator - Clamshell style

LS-9030
Allen® blue boots

LS-9035
Steris/Amsco grey boots

LS-9038
Steris Powerlift - Clamshell style

LS-9040
Allen® Pal-Pro / Foot-Hugger™ style

LS-9045
Great White Legholders pads
Clamshell style

LS-9080
Steris Bariatric Legholder pads
Clamshell style

Gel legholder pads (see pg. 18)
Disposable Legholder pads (see pg. 38)

Extension Pads

LS-9501-2
LS-9501 leg extension pad

LS-9590-2
Amsco leg extension pad

LS-9620-2
Shampaine Radiographic Leg Extension pad 27"L.

For all other makes and models of patient transfer boards and I.A. Extension boards, please call customer service.

Orthopedic Traction Boot Pads

For Orthographic I, Orthographic II.
(sold as pairs only)

LS-9950
For large Tenzor boots

LS-9955
For bariatric Tenzor boots

LS-9980
Secures foot with three deluxe straps.
Attaches to foot plate of traction unit assembly.
Made of soft vinyl with foam padding.

LS-9990
OSI Traction Boot Set - Small

LS-9992
OSI Traction Boot Set - Large

LS-9995
OSI Traction Boot Set - Bariatric X-Large

Knee Crutch Pads

Replacement pads available for most makes and models of knee holders or knee crutches.
Call customer service for more info.

Legholder Pads



LS-9010



LS-9020



LS-9020-D



LS-9025



LS-9028



LS-9030



LS-9038



LS-9040



LS-9080

Orthopedic Traction Boots



LS-9980



LS-9950
LS-9955



LS-9990
LS-9992
LS-9995

Multiple Sizes & Styles available; High quality viscoelastic polymer gel. Provides Anti-shear protection and ultimate comfort in patient positioning.

Head Donuts/Pads

- GP-1500**
Adult (8" Dia.)
- GP-1580**
X Large (8" Dia. x 3")
- GP-1600**
Child (6" Dia.)
- GP-1650**
Infant (3-1/4" Dia.)
- GP-1700**
Adult Horseshoe (8" Dia.)
- GP-1750**
Child Horseshoe (6" Dia.)
- GP-1760**
Infant horseshoe (3.5" Dia.)
- GP-1780**
X large horseshoe (8" Dia. x 3")

Head & Foot Pads

- GP-1800**
Head pillow w/ centering dish (9" x 10" x 1")
- GP-1850**
Head pillow w/ centering dish (9" x 10" x 2")
- GP-1900**
Contoured headrest (7.5" x 7.5" x 2.5" thk)
- GP-1950**
Head pad w/ centering dish (9" x 11" x 3")
- GP-2200**
Anesthesia Prone Headrest w/ tube channels
- GP-2250**
Anesthesia Prone headrest w/ tube channels (9" x 7.5" x 5")
- GP-2800**
2 pc Contoured Ankle gel pads (4" x 7.5" x 3")
- GP-2850**
2 pc Heel gel pads - (3" x 6" x 2.25") w/ molded heel cup indents

Arm Support Pads

- GP-2120**
Armguard Pad (10" x 16")
- GP-2300**
Armboard Pad (6" x 15")
- GP-2325**
Armboard Pad (6" x 20")
- GP-2350**
Armboard Pad (6" x 26")
- GP-2400**
Armboard Pad (4" x 24")
- GP-2600**
Elbow Nerve Protector w/ velcro

Leg Support Pads

- GP-3020**
Knee crutch pad (9" x 12")
- GP-3120**
Knee crutch pad (10" x 16")
- GP-3200**
Knee Crutch Pad (20" x 20")
- GP-3250**
Universal Stirrup Boot Pad (Pair)
- GP-3300**
Candy Cane Stirrup Post Pad (Pair)
- GP-3350**
Candy Cane Ankle Straps (Pair)

Multi-purpose Positioner

- GP-3500**
(6.5" x 3" x 1.75")
- GP-3520**
(13" x 3" x 2")
- GP-3600**
(10" x 5" x 2")
- GP-3700**
(15" x 5" x 2.25")

- GP-3750**
(21" x 5" x 2.5")

Flat Bottom Chest Roll

- GP-4400**
(12" x 7" x 6")
- GP-4420**
(12" x 4" x 4")
- GP-4460**
(16" x 4" x 4")

- GP-4500**
(16" x 7" x 6")

- GP-4600**
(21" x 7" x 6")

- GP-4620**
(12" x 6" x 4")

- GP-4660**
(16" x 6" x 4")

- GP-4680**
(20" x 6" x 4")

- GP-4800**
(6" x 3" x 3")

- GP-4900**
(6" x 2.5" x 2")

- GP-4920**
(10" x 2" x 2")

Surgical Table Pads

- GP-1000**
3/4 Length Table Pad (20" x 46")
- GP-1300**
3 pc Gel Surgical Table Set 1/2" thick without cutout
- GP-1350**
3 pc Gel Surgical Table Set 3/8" thick with cutout

- GP-1400**
Table / Hip Pad (20" x 20")

Specialty Pads

- GP-6100**
Fracture Table Hip Pad w/ 1.75" hole to accommodate perineal post


- GP-6150**
Fracture Table Post Pad (8.25" x 7" x 1/2")

- GP-6300**
Bassinet Gel Pad Overlay (12" x 26")


- GP-6500**
Wilson Frame Pad overlay set w/Velcro 2 pc

- GP-6900**
Bed Overlay - Gurney/Patient bed overlay 34" x 76" x 3/8" thk


Head Donuts/ Pads




Arm Support Pads



Multi-Purpose Positioner



Surgical Table Pads




Head & Foot Pads



Leg Support Pads



Flat Bottom Chest Roll



Specialty Pads





Bean Bag Patient Positioners

Bean bag positioner with replaceable button release valve, durable stain resistant vinyl. Quick and easy hook up to a suction line to evacuate air and form positioner to desired shape. Open valve for positioner to return to relaxed state.

VP-8002
Replacement Valve for bean bag positioners

VP-8010
10" x 20" Positioner

VP-8010G
10" x 20" Positioner w/gel overlay

VP-8020
20" x 20" Positioner

VP-8020G
20" x 20" Bean bag positioner w/gel overlay

VP-8030
30" x 40" Positioner w/Shoulder cut-out

VP-8030G
30" x 40" Positioner w/Shoulder cut-out and gel layer

VP-8036
36" x 40" Positioner w/Shoulder cut-out

VP-8046
46" x 46" Bariatric Positioner w/Shoulder cut-out

VP-8076
30" x 76" Full length Positioner



Ultralite Bag Positioners

- 1/3 the weight of heavy sand
- Re-usable
- Covered with black vinyl
- Great for patient positioning
- Steady & firm support
- Available with or without grommets
- Custom weights & sizes available upon request

SB-2002
2lb. 8" x 11" Positioner

SB-2003
3lb. 9" x 13" Positioner

SB-2005
5lb. 10" x 15" Positioner

SB-2007
7lb. 11" x 16" Positioner

SB-2010
10lb. 12" x 19" Positioner

SB-2202
2lb. 8" x 11" Positioner w/Grommet

SB-2203
3lb. 9" x 13" Positioner w/Grommet

SB-2205
5lb. 10" x 15" Positioner w/Grommet

SB-2207
7lb. 11" x 16" Positioner w/Grommet

SB-2210
10lb. 12" x 19" Positioner w/Grommet



Heavy Fill Sand Bag Positioners

- Re-usable
- Covered with black vinyl
- Great for patient positioning
- Steady & firm support
- Available with or without grommets
- Custom weights & sizes available upon request

SB-1002
2lb. 6.5" x 8" Positioner

SB-1003
3lb. 6.5" x 8" Positioner

SB-1005
5lb. 8" x 10" Positioner

SB-1007
7lb. 8" x 10" Positioner

SB-1010
10lb. 13" x 10.5" Positioner

SB-1202
2lb. 6.5" x 8" Positioner w/Grommet

SB-1203
3lb. 6.5" x 8" Positioner w/Grommet

SB-1205
5lb. 8" x 10" Positioner w/Grommet

SB-1207
7lb. 8" x 10" Positioner w/Grommet

SB-1210
10lb. 10.5" x 13" Positioner w/Grommet

Surgical Table Safety Straps

Safety straps are available in many styles for use on Surgical and other tables.

Rubber Safety Straps w/1 Buckle
Synthetic rubber; Latex free, Durable, and Easy to clean. 2" strap with 1 buckle.

PR-3096
Overall length 96"

PR-3108
Overall length 108"

PR-3120
Overall length 120"

Rubber Safety Straps w/2 Buckles
Synthetic rubber; Latex free, Durable, and Easy to clean. 2" strap with 2 airplane buckles 5"x 20" center panel for patient comfort.

PR-3200
Overall length 96", secures to 56"
Mounts Directly to siderails.

PR-3250
Overall length 96", secures to 56"
Set comes with hooks (PR-0050)



PR-9000



PR-3300
X-tra long, overall length 120"
Secures to 84"

PR-3350
X-tra long, overall length 120"
with hooks (PR-0050)

Safety Strap Mounting Hooks
2 styles to choose from; sold in pairs
Attach to any N. American Surgical table siderail.

PR-0050
Stainless steel hooks, secure onto siderail.

PR-0055
Stainless steel flat hook clips.
Clip on siderail.

Deluxe Safety Straps
4" wide surgical table safety straps,
2 pc Vinyl / Velcro strap; Reusable,
No hooks necessary. Straps secure
directly to siderails.

PR-3400 85"L, overall

PR-3800 104"L, overall

Continuous Safety Straps
4" wide Vinyl with Velcro® closures,
1 pc continuous wrap safety straps. Wraps
completely around table and secure back
on themselves. Ideal for Surgical tables,
Pain Management tables, IA tables, and
Treatment tables.



PR-0055



PR-0050

PR-3920
80" overall length

PR-3950
92" overall length

PR-3980
120" overall length

Alero Armboard Safety Strap PR-4000

2 part strap w/ double Velcro closures
Upper strap 26" L x 2" W. Lower strap
secures around 6" board. Upper strap
secures patient's arm to board.

Economy Armboard Safety Strap PR-4200

6 straps per set; 26" L x 2" W;
1 pc Velcro® strap w/ buckle
Loops around armboard to comfortably
secure patient's arm to armboard.

Stretcher Patient Safety Straps PR-9000

3" wide dura-cloth straps with metal
mounting rings and Velcro closure.
Three 2-pc straps per set.
Straps secure to side of stretcher.

PR-9200
Three 2-pc Straps per set
3" wide, all Velcro strap.



PR-3096



PR-4000



PR-4200





Patient Care Mattresses are designed to work with the patient as well as for the patient. The layers make a conforming surface to the patient's body as well as providing a firm, yet comfortable support. These mattresses are constructed of a combination of solid polyurethane foam and memory foam. And covered in an ultra soft conforming vinyl.

For more details, visit: www.promise-care.com

True Mattress



Model	Dimensions	Type
HBM-1176	34"W x 76"L x 6"T	Standard
HBM-1176RB	34"W x 76"L x 6"T	Roll Barrier
HBM-1180	34"W x 80"L x 6"T	Standard
HBM-1180RB	34"W x 80"L x 6"T	Roll Barrier
HBM-1184	34"W x 84"L x 6"T	Standard
HBM-1184RB	34"W x 84"L x 6"T	Roll Barrier
HBM-1276	36"W x 76"L x 6"T	Standard
HBM-1276RB	36"W x 76"L x 6"T	Roll Barrier
HBM-1280	36"W x 80"L x 6"T	Standard
HBM-1280RB	36"W x 80"L x 6"T	Roll Barrier
HBM-1284	36"W x 84"L x 6"T	Standard
HBM-1284RB	36"W x 84"L x 6"T	Roll Barrier
HBM-1376	38"W x 76"L x 6"T	Standard
HBM-1380	38"W x 80"L x 6"T	Standard
HBM-1384	38"W x 84"L x 6"T	Standard
HBM-1480	40"W x 80"L x 6"T	Standard
HBM-1484	40"W x 84"L x 6"T	Standard
HBM-1580	42"W x 80"L x 6"T	Standard
HBM-1584	42"W x 84"L x 6"T	Standard

True Mattress



- Ultra-Soft 4-way Stretch Vinyl
- Fire Barrier
- Solid Foam Core
- Roll Barrier Panels (Optional upgrade)

Serene Mattress



- Ultra-Soft 4-way Stretch Vinyl
- Fire Barrier
- Solid Foam Core
- Soft Touch Memory Foam
- Dense Support Foam
- Roll Barrier Panel (Optional upgrade)



Serene Mattress

Model	Dimensions	Type
HBM-3176	34"W x 76"L x 6"T	Standard
HBM-3176RB	34"W x 76"L x 6"T	Roll Barrier
HBM-3180	34"W x 80"L x 6"T	Standard
HBM-3180RB	34"W x 80"L x 6"T	Roll Barrier
HBM-3184	34"W x 84"L x 6"T	Standard
HBM-3184RB	34"W x 84"L x 6"T	Roll Barrier
HBM-3276	36"W x 76"L x 6"T	Standard
HBM-3276RB	36"W x 76"L x 6"T	Roll Barrier
HBM-3280	36"W x 80"L x 6"T	Standard
HBM-3280RB	36"W x 80"L x 6"T	Roll Barrier
HBM-3284	36"W x 84"L x 6"T	Standard
HBM-3284RB	36"W x 84"L x 6"T	Roll Barrier
HBM-3376	38"W x 76"L x 6"T	Standard
HBM-3380	38"W x 80"L x 6"T	Standard
HBM-3384	38"W x 84"L x 6"T	Standard
HBM-3480	40"W x 80"L x 6"T	Standard
HBM-3484	40"W x 84"L x 6"T	Standard
HBM-3580	42"W x 80"L x 6"T	Standard
HBM-3584	42"W x 84"L x 6"T	Standard

Serene3 Mattress

Model	Dimensions	Type
HBM-2176	34"W x 76"L x 6"T	Standard
HBM-2176RB	34"W x 76"L x 6"T	Roll Barrier
HBM-2180	34"W x 80"L x 6"T	Standard
HBM-2180RB	34"W x 80"L x 6"T	Roll Barrier
HBM-2184	34"W x 84"L x 6"T	Standard
HBM-2184RB	34"W x 84"L x 6"T	Roll Barrier
HBM-2276	36"W x 76"L x 6"T	Standard
HBM-2276RB	36"W x 76"L x 6"T	Roll Barrier
HBM-2280	36"W x 80"L x 6"T	Standard
HBM-2280RB	36"W x 80"L x 6"T	Roll Barrier
HBM-2284	36"W x 84"L x 6"T	Standard
HBM-2284RB	36"W x 84"L x 6"T	Roll Barrier
HBM-2376	38"W x 76"L x 6"T	Standard
HBM-2380	38"W x 80"L x 6"T	Standard
HBM-2384	38"W x 84"L x 6"T	Standard
HBM-2480	40"W x 80"L x 6"T	Standard
HBM-2484	40"W x 84"L x 6"T	Standard
HBM-2580	42"W x 80"L x 6"T	Standard
HBM-2584	42"W x 84"L x 6"T	Standard

Serene3 Mattress



- Ultra-Soft 4-way Stretch Vinyl
- Fire Barrier
- Solid Foam Core
- Soft Touch Memory Foam
- Inner Core Support Foam Channel
- Flex Base Foam
- Roll Barrier Panel (Optional upgrade)



LS-3000-SR



LS-3500



Lift Assist Leg Positioning System

LS-3000-SR

Yellofins; Completely remanufactured including total tear down, all surfaces refinished, polished, and all wearing parts repaired or replaced.

System Includes:

- 1 set of Pneumatic lift assist leg holders
- 1 set of Pressure management boot pads
- 1 set of siderail clamps

LS-3500

Pneumatic lift assist System
120° degree rotation
Handle controlled abduction and adduction
Mimics natural hip motion
Twist handle to position leg
500 lbs. weight capacity
Gel boot pads available

System Includes:

- 1 set of Pneumatic lift assist leg holders
- 1 set of Pressure management boot pads
- 1 set of siderail clamps

LS-3500-RE

Completely remanufactured including total tear down, all surfaces refinished, polished, and all wearing parts repaired or replaced.

System Includes:

- 1 set of Pneumatic lift assist leg holders
- 1 set of Pressure management boot pads
- 1 set of siderail clamps

Altima Leg Positioning System

ALS-6200

Solid stainless steel adjustable extension bar provides ultimate positioning possibilities. Altima boot clamps are designed for full range adjusting and multi-positioning.

Super-clamps or Clamp-Tite™ clark sockets recommended for mounting legholders to surgical table (sold separately, see pg 42)
Replacement pads & straps available
Gel boot pads available

System Includes:

- 1 pair of Altima Legholders
- 1 set of PM pads Clamshell Style; (ALS-6025)
- 1 set of disposable pads

ALS-6020

Pressure Management Pads (pair)
Limited stock.

ALS-6050

Reusable Gel Pads (Pair)

ALS-6070

Disposable Foam Pads
(10 Pairs per Case)

Universal Well-Leg Support

ALS-6610

Solid stainless steel adjustable extension bar. Boot clamp is designed for full-range positioning. Ideal for positioning patients of any height. Requires one Super clamp or Clamp-Tite™ clark socket for mounting to surgical table (sold separately, see pg 42)

System Includes:

- 1 well Leg Support
- 1 Pressure Management pad

Total Knee Stabilizer

LS-2800

Total Knee Stabilizer, Radiolucent Leg support, Stainless steel, adj. 16" flat rod w/ 10" L support bar covered with 3.5" dia foam pad.
(Replacement pads available. See pg 38)

Knee Crutch Supports

LS-1700

Height adjustable; 360° full rotation
Stainless steel support bars
Includes pressure management pads and safety straps.
Attaches to surgical table using Clamp-tite clark sockets (sold separately, see pg 42)

Complete System Includes:

- 1 pair of Legholders
- 1 set of Safety Straps
- 1 set of Pressure Management pads

Arthroscopic Stress Post

LS-9050-12 - 12" Height Adj.

LS-9050-18 - 18" Height Adj.

Provides a fulcrum to leverage the leg to gain better access to the knee. 5/8" dia. height adjustable post, 4" x 6" pad included. Stainless steel plate & support bar. Available in 12" or 18" post height. Securely attaches to any O.R. Table using a Clamp-Tite™ clark socket (sold separately, see pg 42)



ALS-6050



ALS-6070

Altima Leg Positioning System



ALS-6200
Shown with ALS-6025

Universal Well-Leg Support



ALS-6610

Knee Crutch Supports



LS-1700

Total Knee Stabilizer



LS-2800

Arthroscopic Stress Post



LS-9050-18

LS-9050-12

Candy Cane Style Lithotomy Legholder

LS-2000

Sold in pairs.
Fits all surgical tables.
Adjustable height of 26" - 39"
Support extension 'cane' 10" W.
Stainless steel 'cane' and mounting rod.
Superior flexibility for Lithotomy positioning.
Clamp-Tite™ clark sockets recommended for mounting legholders to surgical table.
(Sold separately See page 42)
Set includes:
1 pair of ankle support straps

Lithotomy Legholder Ankle Straps

LS-2020

Durable cotton straps
Double stitched for superior strength.
Sold in pairs.

Superior Strength Lithotomy Legholder Heel Cups

LS-2040

Extra soft heel support cups
Sold in pairs.

Lithotomy Legholder Gel Post Pads

GP-3300

Viscoelastic Gel pads slip over lithotomy post to protect patient extremity.
Sold in pairs.

Lithotomy Legholder Gel Ankle Straps

GP-3350

Comfortable support
Sold in pairs.

Arthroscopic Leg Holder Unit

LS-1400

Accommodates all patients.
Adjusts medially and laterally to apply desired pressure.
Dual hinged positioning clamp.
Spring loaded adjusting lever.
Siderail clamp and pad included.

Replacement Pads for Arthroscopic Legholder

LS-1420

2-ply foam pad
10 pads per case

Legholder/ Accessory Cart

LS-1180

Solid stainless steel Legholder Cart
Ideal for storage of legholder system and other O.R. accessories.
Casters allow easy mobility.



LS-1180

Pickett Fence Leg Prep Unit

LS-4200

Efficiently preps both legs
5/8" dia. Stainless steel support bar
Height adjustable up to 30"
Attaches to surgical table using a Clamp-Tite clark socket (Sold separately See page 42)

LS-4300

Bariatric Pickett Fence Leg Prep unit
Easy positioning
Efficiently preps both legs
5/8" dia. Stainless steel support bar
Height adjustable up to 30"
Attaches to surgical table using a Clamp-Tite clark socket (Sold separately. See page 42)

Uni-Prepper Leg Prep Unit

LS-4800

Most versatile leg prepper available
Infinite positioning possibilities
Leg cradle measures 8"W x 5"H x 3.5"D.
Adjustable height of 12" - 27"
Stainless steel construction. Fits all surgical tables. Pressure management pad included.
Requires a Clamp-Tite clark socket for mounting (Sold separately. See page 42)

Bariatric Uni-Prepper Leg Prep Unit

LS-4850

Equivalent to LS-4800 with heavier construction and 10.5" wide support cradle.
Requires a Clamp-Tite clark socket for mounting (Sold separately. See page 42)

Replacement Uni-Prepper Pad

LS-4801-1

Memory foam with 4-way stretch vinyl pad.
For LS-4801

Replacement Bariatric Uni-Prepper Pad

LS-4851-1

Memory foam with 4-way stretch vinyl pad.
For LS-4851

Candy Cane Style Lithotomy Legholder



LS-2000

LS-2020

LS-2040

GP-3300

GP-3350


Arthroscopic Leg Holder Unit



LS-1400

LS-1420

Pickett Fence Leg Prep Unit



LS-4300

LS-4200

Uni-Prepper Leg Prep Unit



LS-4800

Bariatric Uni-Prepper Leg Prep Unit



LS-4850

Leg/Foot Extension

LS-9500

Stainless steel and phenolic construction. Adds ten inches to length of O.R. table. Doubles as a foot support in vertical position. Comes with 2" pad and mounting clamps. Quick and easy to setup. This extension is designed for most surgical tables – Specify make and model of surgical table when ordering.

Stainless Steel Straight Bar Accessory Clamp

CL-1000

Ergonomic handle. Solid stainless steel construction. Accommodates all rectangle bar accessories up to 1/4" x 1". Easily attaches to siderails of all North American surgical tables.

CL-1100

Solid stainless steel construction. Accommodates all rectangle bar accessories up to 1/4" x 1". Easily attaches to siderails of all North American surgical tables.

CL-1500

Aluminium construction. Accommodates all rectangle bar accessories up to 1/4" x 1". Easily attaches to siderails of all North American surgical tables.



LS-9500
Shown without pad

Clamp-Tite™ Stainless Steel Clark Socket

CL-8000-SS

Design ensures ultimate clamping ability for all surgical table accessory positioning needs. Accommodates all accessories up to 3/4" Dia. Attaches to siderails of all North American surgical tables. Allows secure positioning of surgical table accessories.

Standard Clark Socket

CL-8100

Standard Style Clark Socket. Solid Stainless Steel Construction. Accommodates round bar accessories up to 3/4" dia. Attaches to all North American surgical table side rails.

Stainless Steel Round Bar Accessory Clamp

CL-2000

Ergonomic Handle. Solid stainless steel construction. Accommodates all round bar accessories 1/4" up to 3/4" Dia. Easily attaches to siderails of all North American surgical tables.

CL-2100

Solid stainless steel construction. Accommodates all round bar accessories 1/4" up to 3/4" Dia. Easily attaches to siderails of all North American surgical tables.

Super Clamp Multi-Positioner

CL-7500

Solid stainless steel construction. Effortlessly secures most accessories to siderail of any North American surgical table. Both round bars up to 3/4" Dia. as well as straight bars up to 1/4" x 1". Also accommodates square bar accessories up to 3/4" sq.

Optional: Super Clamp to accommodate European siderails available (CL-7500-E)

Siderail Extender

EX-1000

Extends the siderail of surgical tables by 8" to accommodate more mounting accessories. Provides great patient positioning flexibility Solid Stainless Steel. Doubles as a straight bar clamp. Accommodating all straight bar accessories up to 1" x 1/4". Attaches to all North American surgical table siderails. Sold in pairs.

Table Width Extender

EX-2200

Surgical table side/width extension. Easily mounts to existing table siderail. Extends table surface 4" x 20". Comes with 2" pad and 18" Accessory Siderail.



CL-7500

Square Bar



Round Bar



Flat Bar

Stainless Steel Straight Bar Accessory Clamp

CL-1000



CL-1100



CL-1500



Clamp-Tite™ Stainless Steel Clark Socket

CL-8000-SS



Standard Clark Socket

CL-8100



Stainless Steel Round Bar Accessory Clamp

CL-2000



CL-2100



Siderail Extender



EX-1000

Table Width Extender



EX-2200

Patient Monitors

DATASCOPE PASSPORT V

The Passport V offers a comprehensive range of options that make it well-suited for most hospital departments. The flexible feature-set includes your choice of CO₂ and SpO₂ technologies. [See website for more product details...](#)

DATEX OHMEDA CARDIOCAP 5

- Small footprint bedside monitor with large 10.4" active matrix color display
 - Proven easy-to-use S/5 user interface
 - Shows up to 6 waveforms
 - Color-coded, prioritized alarms
 - 4 user-configurable modes for care-specific preconfigurations
- [See website for more product details...](#)

DINAMAP MPS SELECT

- With a variety of configurations, the Dinamap MPS Select can help you meet the needs of your medical facility.
 - It is compact for maximum portability.
 - Available parameters include NIBP, SpO₂, ECG and Temp.
 - Weighs less than 10 lbs.
 - Optional printer available
 - AC or battery power operation
 - Audible low battery indicator to prevent memory loss
 - Optional spot check SpO₂
 - Reliability and performance in one Easy to use interface
- [See website for more product details...](#)

GE CARESCAPE B650

- Integrates with GE Engstrom and Anesthesia Carestation devices with remote on/off function.
 - Monitoring can start from the anesthesia station.
 - CARESCAPE Patient Data Module offers seamless patient data capture in HIS/ EMR for continuity of care.
 - Key parameters help with quality of care: dual SpO₂, continuous metabolics and gas exchange, Ek-Pro™ arrhythmia analysis, and Adequacy of Anesthesia parameters, such as Entropy and NMT.
- [See website for more product details...](#)

GE CARESCAPE B850

The GE CARESCAPE B850 Patient Monitor is a high-acuity monitor with various parameters to help clinicians assess the adequacy of anesthesia in the operating room. Manage patient flow efficiently with the monitor's data continuity and integration across care areas. [See website for more product details...](#)

GE DASH 3000 PRO

- Dash 3000 is specially built to withstand the punishment of your most demanding intra-hospital transport applications with certifications for both ground and air transport
 - Compact: Weighing no more than 12 pounds and measuring only 8½" x 10½" x 8" regardless of configuration, these compact monitors make transport quick and trouble free
- [See website for more product details...](#)

GE DASH 5000

- Advanced clinical parameters including 12SL™ for simultaneous ECG monitoring, GE EK-Pro™ for multi-lead arrhythmia analysis, ICG, mainstream and sidestream CO₂ and GE DINAMAP® NIBP with SuperSTAT™ Networking options include wired and wireless
 - High resolution CRG trends for NICU monitoring
- [See website for more product details...](#)

GE PROCARE B40 PATIENT MONITOR

The GE Procure B40 monitor utilizes advanced GE algorithms and proven clinical measurement technologies for accurate and reliable monitoring even in challenging conditions. The B40 is suitable for monitoring adult, pediatric and neonatal patients. It is easily configured to monitor a variety of possible parameters including NIBP, IBP, ETCO₂, heart rate, respiration rate, ECG, temperature, SpO₂, agent identification, in addition to gas and respiratory parameters. [See website for more product details...](#)

PHILIPS INTELLIVUE MP20

The IntelliVue MP20 patient monitor provides powerful monitoring capability and essential measurements in a compact package to match the pace and unique needs of immediate care, special procedures, lower acuity environments, and patient transfer. The user interface has been redesigned to improve visibility of patient data, making it easier to use, and to enhance compatibility with standard software. [See website for more product details...](#)

PHILIPS INTELLIVUE MP30


The IntelliVue family of networked patient monitors gives care teams throughout the hospital more of the information they need right at the patient's side. The Philips IntelliVue MP30 provides powerful monitoring capability and essential measurements in a compact package to match the pace and unique needs of intermediate care, ambulatory surgery, post-operative care, special procedures, lower acuity environments, and patient transfer. Featuring highly flexible screen configurations designed to suit patient acuity, department protocols, or specific procedure requirements. [See website for more product details...](#)

PHILIPS INTELLIVUE MP5

Philips IntelliVue MP5 patient monitor provides a clear connection to your patients, delivering IntelliVue monitoring power and functionality in a compact, rugged housing to match the demands of a wide range of care environments in and out of the hospital.

The IntelliVue MP5 offers advanced features such as an intuitive touchscreen with clinical measurements, crisp and clear display, and one-touch commands. Clinical Decision Support tools including Horizon Trends, ST Map, and ProtocolWatch reveal more information when it counts most, enhancing clinician confidence. [See website for more product details...](#)

Patient Monitors



DATASCOPE PASSPORT V



DATEX OHMEDA CARDIOCAP 5



DINAMAP MPS SELECT



GE CARESCAPE B650




GE CARESCAPE B850



GE DASH 3000 PRO



GE DASH 5000



GE PROCARE B40 PATIENT MONITOR



PHILIPS INTELLIVUE MP20



PHILIPS INTELLIVUE MP30



PHILIPS INTELLIVUE MP5

Surgical Tables

BERCHTOLD OPERON® OR TABLE D850, D820, D760

The Bертold OPERON® Modular OR Table comes in three varying bases, the D 850, the D 820, and the D 760. Each base offers varying height adjustment and positioning surfaces.

There are a wide variety of accessories available for the Bертold OPERON® to accommodate multiple surgical modalities.

[See website for more product details...](#)

D 850

Features

- Height Adjustment: 22 - 46 in (580 - 1180 mm)
- Maximum Patient Weight: 1250 lbs (567 kg)
- Longitudinal Slide Ideal for seated or standing surgical procedures
- Maximum flexibility for imaging
- Foot pump safety system and backup hand pendant

D 820

Features

- Height Adjustment: 25 - 42 in (635 - 1075 mm)
- Maximum Patient Weight: 1000 lbs (454 kg)
- Longitudinal Slide
- Ideal for multiple inpatient and outpatient surgeries
- Foot pump safety system and emergency hand pendant under table
- Maximum flexibility for imaging

D 760

Features

- Height Adjustment: 25 - 42 in (635 - 1075 mm)
- Maximum Patient Weight: 1000 lbs(454 kg)
- No Longitudinal Slide

SKYTRON 6500 HD O.R. TABLE

Skytron 6500 HD surgical tables are available in manual, electric, and remote control models. In addition to the 6500 HD, refurbished models are available from the Skytron 5000, 6000, 6500, and 7000 series.

Features

- General use electro-hydraulic surgical table
- Can lift and fully articulate patients weighing up to 850 lbs.
- C-Arm compatible
- X-ray top set
- 180 degree table top rotation
- Return to level
- Buttons are clearly labeled, making the table easy to use and enhancing patient safety
- Fail-safe feature makes sure brakes are locked prior to table use
- Stainless steel cover is attractive and easy to clean
- Self-leveling, four-point brake system
- Built-in kidney bridge
- Manually-operated head section features a full width quick-release bar
- Removable powered leg section
- Narrow base
- Positions include:
 - Trendelenburg: Back (raise or lower); Side (lateral tilt); Foot (raise or lower); Flex ("v" or inverted "v")

[See website for more product details...](#)

SKYTRON 6500 O.R. TABLE

The Skytron 6500 Elite General Purpose Surgery Table includes operating hand control and new internal battery. It is in excellent condition and certified to proper operating specifications. Warranty: 1 Year - Parts and Labor. The 6500 further extends the range of applications of the 6000 series. It increases the patient weight limit to 500 lb, as well as increases the degrees of articulation offered. It also incorporates a powered kidney bridge as standard, along with fully

powered brakes and a "return-to-level" function, all of which are controlled from the hand control. The rotating top of 210° further enhances C-arm and microscope access, along with more leg room for a seated surgeon and in many cases, better surgical access. A stable footprint and rigid surface provide assurance that the patient position will be maintained as desired. The 6500 model features a detachable head section that can be inserted into the end of the foot section to provide support for legs when positioning patient for lithotomy procedures. It is also equipped with a manual brake release that is independent of the main controls and requires no power to operate.

[See website for more product details...](#)

SKYTRON 6701 SURGERY TABLE

The Skytron Hercules 6701 is a surgery table that can be used for general surgery or for bariatric patients.

Features

- 1,000-lb. articulate capacity
- 1,200 lb. lift
- Auto beach chair position
- 41-inch maximum height
- 23-inch minimum height
- Power kidney bridge
- 41-inch maximum height
- Removable back and leg section

[See website for more product details...](#)

Surgical Tables



BERCHTOLD OPERON® OR TABLE - D850, D820, D760



SKYTRON 6500 HD O.R. TABLE



SKYTRON 6500 O.R. TABLE



SKYTRON 6701 SURGERY TABLE



800-277-5852
WWW.AMGEQUIPMENT.COM