





Engineered and Manufactured in the U.S.

## **Sure-Dry Inspiration**

The Sure-Dry Scope Cabinet includes the industry best Positive Pressure HEPA Filtered Drying system. Inside air is replaced inside the cabinet 4 times per minute. Meets or exceeds Joint Commission, AORN, SGNA, and AAMI guidelines.



## SCSS8136DRDP (shown above) Stainless Steel 16 Scope Capacity with Key Lock. Cabinet also available in powder coat painted low carbon steel (16 colors available).

## **Scope Cabinet Features:**

- · Fully redundant fan system
- Industries only unit with 3 intake and 2 exhaust fans
- Fan run time is fully adjustable and activated each time cabinet is closed
- HEPA filter removes 99.7% of particulates, .03 Microns
- Holds 16 scopes with secure adjustable mounting system
- Replaceable HEPA Filter system without removal of scopes
- Medical grade electrical system
- Tempered glass double doors
- Multiple locking options available including badge entry
- Removable drip tray
- Cabinets are delivered fully assembled
- Harloff's Industry Best 12 year limited manufacturers warranty
- Overall cabinet dimensions: 93"H x 36"W x 21"D

Simple "Plug and Play" upgradeable with Dri-Scope Aid®, industry leading internal lumen drying system