





Engineered and Manufactured in the U.S.

Sure-Dry Inspiration

The Sure-Dry Scope Cabinet includes the industry best Positive Pressure HEPA Filtered Drying system. Inside air is replaced inside the cabinet 4 times per minute. Meets or exceeds Joint Commission, AORN, SGNA, and AAMI guidelines.



SCSS8136DREDP (shown above)

Stainless Steel 16 Scope Capacity with basic electronic lock. Cabinet also available in powder coat painted low carbon steel (16 colors available).

Scope Cabinet Features:

- Fully redundant fan system
- Industries only unit with 3 intake and 2 exhaust fans
- Fan run time is fully adjustable and activated each time cabinet is closed
- HEPA filter removes 99.7% of particulates, .03 Microns
- Holds 16 scopes with secure adjustable mounting system
- Replaceable HEPA Filter system without removal of scopes
- Medical grade electrical system
- Tempered glass double doors
- Multiple locking options available including badge entry
- Removable drip tray
- Cabinets are delivered fully assembled
- Harloff's Industry Best 12 year limited manufacturers warranty
- Overall cabinet dimensions: 93"H x 36"W x 21"D

Simple "Plug and Play" upgradeable with Dri-Scope Aid®, industry leading internal lumen drying system