

# **Generation 3 Galaxy 4X4**

Mobile Exam Light

Part #: 4X4-GSM

### Introduction

The **StarTrol** Galaxy Generation 3, the newest member of our lighting family. The LED's are brighter than our current lights and features five level diming. The new sleek design, features a fully enclosed light head for improved clean ability. **StarTrol** is proudly made in the United States, with an unmatched 10 year limited warranty.

#### **Standard Features**

Touch Sense on / off with 5 step dimming, along with level indicator lights

Multiple configurations: Mobile, Wall, Single, and Dual Ceilings feature spring arm suspension Removable anodized aluminum handle Enclosed light head for improved clean ability 10 year limited warranty

# **Architectural Specifications**

**StarTrol** Gen 3 4X4-GSM has an 8.75" diameter aluminum light head featuring 16 LEDs (58,000 lux at a meter). This light has an arm reach of 31". The light head has a 700° rotation at the arm and head pivot point and 330° at the light head joint. The overall height of the light is 60.5", the dimensions of mobile base are 21" x 20.1" x 12.5". Electrical power is supplied by a 10' hospital grade power cord (included).

# **Technical Data**

Part number: 4X4-GSM

Input voltage: 100-240 VAC 50-60Hz Auto Adjusting

**Input Current:** .33 Amps @ 120 VAC **Light Output:** 58,000 lux at a meter

**Spot Size:** 8 inches

**Color Temperature:** 4500K **Color rendering index:** 95

**R9:** 99

**LED Life:** 50,000+ hours **Number of LEDs:** 16

**Dimming:** 5 Level Touch Sense

**Total Reach:** 31 inches

**Base dimension:** 21" X 20.1" 12.5"



# **Accessories**

Additional Aluminum handle (Part# HDN) Replacement Power cord (Part# PWC)