

Generation 3 Galaxy 8X4

Dual Ceiling Mounted Light

Part #: 8X4-GS2



Introduction

The **StarTrol** Galaxy Generation 3, the newest member of our lighting family. The LED's are brighter than our current lights and features five level diming. The new sleek design, features a fully enclosed light head for improved clean ability. **StarTrol** is proudly made in the United States, with an unmatched 10 year limited warranty.

Standard Features

Unlimited Rotation at the 2 main pivot points Robust spring - that never needs adjustment

Touch Sense on / off with 5 step dimming, along with

level indicator lights

Multiple configurations (Wall, Single, and Dual Ceilings)

Removable anodized aluminum handle

Enclosed light head for improved clean ability 10 year limited warranty

Accessories

Additional Aluminum handle (Part# HDN)
Ceiling Cover for flush mount installations
(Part # GHC)

Jumbo Bezel for above ceiling installations (Part# JZB)

Architectural Specifications

StarTrol Gen 3 8X4-GS2 is made up of two 12.5" diameter aluminum light heads featuring 32 LEDs per head (200,000 lux at a meter). This light has an arm reach of 120" with unlimited rotation at the 2 main pivot points. The light head has a 700° rotation at the arm and head pivot point and 330° at the light head joint.

Technical Data

Part number: 8X4-GS2

Input voltage: 100-240 VAC 50-60Hz Auto Adjusting

Input Current: 2 Amps @ 120 VAC **Light Output:** 200,000 lux at a meter

Spot Size: 10 inches

Color Temperature: 4500K **Color rendering index:** 95

R9: 99

LED Life: 50,000+ hours **Number of LEDs:** 64

Dimming: 5 Level Touch Sense **Total Reach:** 60 inches (5 feet)

Minimum Ceiling Height: 96" with an 80" pivot** **Maximum Ceiling Height:** 126" with an 80" pivot**

^{**} Check with the manufacturer for lights outside these measurements.