# 3-Drawer Junior ShortIsolation Cart (JSGKU-399-YEL)

👿 waterloohealthcare.com/product/3-drawer-junior-short-jsgku-399-yel



Economical steel isolation carts offer a robust work platform that can be configured with accessories to address your code cart needs. Shares many of the same great features as our aluminum cart but at a more economical initial price. Other features include:


Smooth ball bearing slides with auto-close mechanism to pull the drawer fully closed.

Advanced nylon reinforced polycarbonate bumper provides a wider wheel base and greater cart stability and protection. Single piece construction resists scuffs and scratches.

5" Premium "Easy roll" sealed bearing casters with precision bearing swivels. Thread guards minimize sutures and other materials from collecting around axles and keep carts rolling smoothly.

# Description

# SPECIFICATIONS



## Frame Color: LTG



### Drawer Color: YELLOW Light Gray LTG Natural NTL White WHT Yellow YEL Red RED Wine WNE Pink PNK Orange ORG Dark Blue DKB Electric Blue ELB Violet VIL Light Blue LTB Teal Green TLG Spring Green SPG Fairway Green FWG Dark Gray GRA Wood Grain Finishes Natural Wood NAT Cherry CHY Walnut WNT


# Key lock Lever lock Image: Set lock Push button lock Image: Set lock I

Lock System: KEY

LOCK

Dual Access lock

# Top: REMOVEABLE PLASTIC



Flexible plastic top

# Bumpers: STABILIZER BASE







# Casters: 5" STANDARD



