

## 3-Drawer Short (ESGKU-399-LTG)

[waterloohealthcare.com/product/3-drawer-short-esgku-399-ltg](https://waterloohealthcare.com/product/3-drawer-short-esgku-399-ltg)



Waterloo value carts are designed with budget and size in mind. If you need a compact mobile storage solution that keeps an eye on your pocketbook then our Value cart is for you!

Standard Features include:

- Centralized key lock

- Ball bearing glide drawers.

- Removable ABS plastic top.

- 5" Premium easy roll casters (two braking, one steering, and one full swivel).

- Stabilizer base with high impact nylon reinforced polycarbonate bumper.

SKU: ESGKU-399-LTG Category: Economy/Value Carts

## Description

---


## SPECIFICATIONS

---

Frame Type: STEEL-JUNIOR  
Frame Size: SHORT

Drawer  
Dimensions: 16.5"w x 9-9  
16.5"d

Drawer Height: 3-



3"  1  
drawer

6" 


9"  2  
drawers

12" 

Dimensions: 23"w x 24.5"d x  
36"h

---

Frame Color: LTG



Wood Grain Finishes



Drawer Color:  
LIGHT GRAY



Lock System: KEY  
LOCK



Top: REMOVEABLE PLASTIC



Flexible plastic top

Bumpers:  
STABILIZER BASE



Stabilizer Base



No bumpers

Casters: 5"  
STANDARD



5" Caster (1 locking)



5" Caster (1 locking)