### 3-Drawer ShortIsolation Cart (USGKU-399-YEL)

**waterloohealthcare.com**/product/3-drawer-short-usgku-399-yel



Economical steel isolation carts offer a robust work platform that can be configured with accessories to address your code cart needs. Shares many of the same great features as our aluminum cart but at a more economical initial price. Other features include:

Smooth ball bearing slides with auto-close mechanism to pull the drawer fully closed.

Advanced nylon reinforced polycarbonate bumper provides a wider wheel base and greater cart stability and protection. Single piece construction resists scuffs and scratches.

5" Premium "Easy roll" sealed bearing casters with precision bearing swivels. Thread guards minimize sutures and other materials from collecting around axles and keep carts rolling smoothly.

SKU: USGKU-399-YEL Category: <u>Isolation Carts</u>

# Description

## **SPECIFICATIONS**

Frame Type: STEEL Frame Size: SHORT

Drawer Dimensions: 22"w x

9-9


Drawer Height: 3-

16.5"d



Dimensions: 29"w x 24.5"d x

37.5"h



#### Frame Color: LTG



#### Drawer Color: Yellow




### Lock System: KEY LOCK



Dual Access lock

#### Top: REMOVEABLE PLASTIC



Flexible plastic top

#### Bumpers: STABILIZER BASE







No bumpers

#### Casters: 5" **STANDARD**







3" Caster (1 locking)