

4-Drawer Junior X-Tall Isolation Cart (JXGKA-3999-YEL)

 waterloohealthcare.com/product/4-drawer-junior-x-tall-jxgka-3999-yel



This cart is built on the tallest of our Aluminum frames and comes standard with a full set of drawer dividers, transport lever lock and removable plastic top for easy cleaning. Other features include:

Dual slide out work surfaces operate independently and virtually double your work surface!

Smooth ball bearing slides with auto-close mechanism to pull the drawer fully closed.

Large stylish/contemporary drawer handles offer a more updated look and enhanced functionality

Adjustable metal dividers (One long and One short) included in all aluminum cart drawers.

Advanced nylon reinforced polycarbonate bumper provides a wider wheel base and greater cart stability and protection. Single piece construction resists scuffs and scratches.

5" Premium "Easy roll" sealed bearing casters with precision bearing swivels. Thread guards minimize sutures and other materials from collecting around axles and keep carts rolling smoothly.

SKU: JXGKA-3999-YEL Category: Isolation Carts




Description

SPECIFICATIONS

Frame Type: ALUM-JUNIOR
Frame Size: X-TALL


Drawer
Dimensions: 16.5"w x
16.5"d

Drawer Height: 3-
9-9-9



Dimensions: 23"w x 24.5"d x
45"h

Frame Color: LTG



Wood Grain Finishes



Drawer Color:
YELLOW



Lock System: KEY
LOCK



Top: REMOVEABLE PLASTIC



Flexible plastic top

Bumpers:
STABILIZER BASE



Stabilizer Base

No bumpers

Casters: 5"
STANDARD



5" Caster (1 locking)

3" Caster (1 locking)