4-Drawer Tall Isolation Cart (UTGKU-3699-YEL)

waterloohealthcare.com/product/4-drawer-tall-utgku-3699-yel



Economical steel isolation carts offer a robust work platform that can be configured with accessories to address your code cart needs. Shares many of the same great features as our aluminum cart but at a more economical initial price. Other features include:

Smooth ball bearing slides with auto-close mechanism to pull the drawer fully closed.

Advanced nylon reinforced polycarbonate bumper provides a wider wheel base and greater cart stability and protection. Single piece construction resists scuffs and scratches.

5" Premium "Easy roll" sealed bearing casters with precision bearing swivels. Thread guards minimize sutures and other materials from collecting around axles and keep carts rolling smoothly.

SKU: UTGKU-3699-YEL Category: Isolation Carts

Description

SPECIFICATIONS

Frame Type:

Drawer

Drawer Height: 3-6-9-9

STEEL

Dimensions:

22"w x 16.5"d Frame Size: TALL





drawer

drawers

Dimensions:

29"w x

24.5"d x

43.5"h

Frame Color: LTG



Drawer Color: Yellow



Lock System: KEY LOCK





Key lock

Lever lock





Gate lock

Push button lock





Keypad lock

Proximity lock



Dual Access lock

Top: Bumpers: Casters: 5" STANDARD

REMOVEABLE STABILIZER PLASTIC BASE











5" Caster (1 locking)

3" Caster (1 locking)