5-Drawer MediumEmergency Cart (UMRLA-33369-RED)

waterloohealthcare.com/product/6-drawer-x-tall-umrla-33369-red



This cart is built on the medium version of our Aluminum frames and comes standard with a full set of drawer dividers, transport lever lock and removable plastic top for easy cleaning. Other features include:

Dual slide out work surfaces operate independently and virtually double your work surface!

Smooth ball bearing slides with auto-close mechanism to pull the drawer fully closed.

Large stylish/contemporary drawer handles offer a more updated look and enhanced functionality

Adjustable metal dividers (One long and One short) included in all aluminum cart drawers.

Advanced nylon reinforced polycarbonate bumper provides a wider wheel base and greater cart stability and protection. Single piece construction resists scuffs and scratches.


5" Premium "Easy roll" sealed bearing casters with precision bearing swivels. Thread guards minimize sutures and other materials from collecting around axles and keep carts rolling smoothly.

SKU: UMRLA-33369-RED Category: Emergency Carts Tag: Emergency Cart Aluminum

Description

SPECIFICATIONS

Frame Type: ALUM-UNICART Drawer Drawer Height: 3-Frame Size: MEDIUM Dimensions: 22"w x 3-3-6-9 16.5"d



Frame Color: RED




Drawer Color:

RED




Lock System: LEVER LOCK



Dual Access lock

Top: REMOVEABLE PLASTIC



Flexible plastic top

Bumpers: STABILIZER BASE





No bumpers

Casters: 5" STANDARD



