

5-Drawer Short Anesthesia Cart (USGKU-33339-TLG)

 waterloohealthcare.com/product/5-drawer-short-usgku-33339-tlg



Economical steel carts offer a robust work platform that can be configured with accessories to address your anesthesia cart needs. Shares many of the same great features as our aluminum cart but at a more economical initial price. Other features include:

- Smooth ball bearing slides with auto-close mechanism to pull the drawer fully closed.

- Advanced nylon reinforced polycarbonate bumper provides a wider wheel base and greater cart stability and protection. Single piece construction resists scuffs and scratches.
- 5" Premium "Easy roll" sealed bearing casters with precision bearing swivels. Thread guards minimize sutures and other materials from collecting around axles and keep carts rolling smoothly.

SKU: USGKU-33339-TLG Category: Anesthesia Carts



Description

SPECIFICATIONS

Frame Type: STEEL
Frame Size: SHORT

Drawer
Dimensions: 22"w x
16.5"d

Drawer Height: 3-
3-3-9



Dimensions: 29"w x 24.5"d x
36"h

Frame Color: LTG



Drawer Color:
TEAL



Lock System: KEY
LOCK



Top: REMOVEABLE PLASTIC



Bumpers:
STABILIZER BASE



Casters: 5"
STANDARD

