5-Drawer Short (ESGKU-33366-LTG)

waterloohealthcare.com/product/5-drawer-short-esgku-33366-ltg-2



Waterloo value carts are designed with budget and size in mind. If you need a compact mobile storage solution that keeps an eye on your pocketbook then our Value cart is for you!

Standard Features include:

Centralized key lock

Ball bearing glide drawers

Removable ABS plastic top

5" Premium easy roll casters (two braking, one steering, and one full swivel)

Stabilizer base with high impact nylon reinforced polycarbonate bumper

SKU: ESGKU-33366-LTG Category: <u>Economy/Value Carts</u>

Description

SPECIFICATIONS

Frame Type: STEEL-JUNIOR Frame Size: SHORT	Drawer Dimensions: 16.5"w x 16.5"d	Drawer Height: 3- 3-3-6-6
	16.5" Economy size	³ drawers
		97
Dimensions: 23"w x 24.5"d x 36"h		127

Frame Color: LTG



Drawer Color: LIGHT GRAY



Lock System: KEY LOCK Lever lock Key lo






Proximity lock



Keypad lock

Dual Access lock

Top: REMOVABLE PLASTIC



Flexible plastic top

Bumpers: STABILIZER BASE





Casters: 5" **STANDARD**





3/3