6-Drawer Tall Anesthesia Cart (UTGSU-333369-DKB)

👿 waterloohealthcare.com/product/6-drawer-tall-utgsu-333369-dkb



Economical steel carts offer a robust work platform that can be configured with accessories to address your anesthesia cart needs. Shares many of the same great features as our aluminum cart but at a more economical initial price. Other features include:


Smooth ball bearing slides with auto-close mechanism to pull the drawer fully closed.

Advanced nylon reinforced polycarbonate bumper provides a wider wheel base and greater cart stability and protection. Single piece construction resists scuffs and scratches.

5" Premium "Easy roll" sealed bearing casters with precision bearing swivels. Thread guards minimize sutures and other materials from collecting around axles and keep carts rolling smoothly.

SKU: UTGSU-333369-DKB Category: Anesthesia Carts

SPECIFICATIONS




Frame Color: LTG



Drawer Color: DARK BLUE




Lock System: PUSH BUTTON LOCK





Dual Access lock

Top: REMOVEABLE PLASTIC



Bumpers: STABILIZER BASE









3" Caster (1 locking)

Flexible plastic top