6-Drawer X-Tall Anesthesia Cart (UXGSU-333669-ELB)

waterloohealthcare.com/product/6-drawer-x-tall-uxgsu-333669-elb



Economical steel carts offer a robust work platform that can be configured with accessories to address your anesthesia cart needs. Shares many of the same great features as our aluminum cart but at a more economical initial price. Other features include:

Smooth ball bearing slides with auto-close mechanism to pull the drawer fully closed.

Advanced nylon reinforced polycarbonate bumper provides a wider wheel base and greater cart stability and protection. Single piece construction resists scuffs and scratches.

5" Premium "Easy roll" sealed bearing casters with precision bearing swivels. Thread guards minimize sutures and other materials from collecting around axles and keep carts rolling smoothly.

SKU: UXGSU-333669-ELB Category: <u>Anesthesia Carts</u>

Description

SPECIFICATIONS

Frame Type: STEEL
Frame Size: X-TALL

Dimensions: 22"w x 3-3-6-6-9

16.5"d

The state of the sta


Frame Color: LTG



Drawer Color: ELECTRIC BLUE



Lock System: PUSH BUTTON LOCK



Dual Access lock

Top: REMOVEABLE PLASTIC



Flexible plastic top

Bumpers: STABILIZER BASE







No bumpers

Casters: 5" STANDARD





