

9-Drawer X-Tall PEDIATRIC Cart (UXGLA-9PEDS)

waterloohealthcare.com/product/9-drawer-x-tallpediatric-cart-uxgla-9ped



This cart is built on the medium version of our Aluminum frames and comes standard with a full set of drawer dividers, transport lever lock and removable plastic top for easy cleaning. Other features include:

Dual slide out work surfaces operate independently and virtually double your work surface!

Smooth ball bearing slides with auto-close mechanism to pull the drawer fully closed.

Large stylish/contemporary drawer handles offer a more updated look and enhanced functionality

Adjustable metal dividers (One long and One short) included in all aluminum cart drawers.

Advanced nylon reinforced polycarbonate bumper provides a wider wheel base and greater cart stability and protection. Single piece construction resists scuffs and scratches.

5" Premium "Easy roll" sealed bearing casters with precision bearing swivels. Thread guards minimize sutures and other materials from collecting around axles and keep carts rolling smoothly

SKU: UXGLA-9PEDS Category: Emergency Carts Tag: Emergency Cart Aluminum



Description

SPECIFICATIONS

Frame Type: ALUM-PED
Frame Size: X-TALL

Drawer Dimensions:
22"w x 16.5"d

Drawer Height: 3-3-3-3-3-3-3-6



3"  8 drawers

4" 

6"  1 drawer

9" 

12" 

Dimensions: 29"w x
24.5"d x 45"h

Frame Color: LTG



Drawer Color: MULTI-COLOR



Lock System: LEVER LOCK



Top: REMOVEABLE PLASTIC



Bumpers: STABILIZER BASE



Casters: 5" STANDARD

